

Computer Ethernet Rack

Pre-configured and customized instrument PC With full support

Advantages:


- PEAK compatibility
- One year remote support* including advice on plug-ins, e.g. manipulator control
- Win10 64-bit computer with hardware and software necessary for analyser control
- 'Plug and Play', customized and pre-installed
- Most recent SES software version with Real-time straight slit compensation
- Igor Pro 8 Software License



Standard PC

DA30 PC

A Pre-configured and customized computer that meets the new PC requirements for PEAK compatibility. Latest platform; drivers and SES installed and tested. Including one year software support. Send instrument files to services-e-spec@Scientaomicron.com



Summary

Part numbers:

- 201042 Computer SES ethernet rack
- 201043 Computer SES DA30 ethernet rack

Parts included:

- PC (does not include monitor/display)
- Power cable
- Instruction

Hardware prerequisites:

- 102133 DOM upgrade
- 70581XX Card HVC20 (with firmware ≥ 3.3)
- Depending on current configuration, additional upgrades may be required

*Prerequisite for included remote support:

- Complete Ethernet upgrade solution, see options

Options:

- 201040 Ethernet Rack upgrade (required for PEAK compatibility and remote support)
- 85827/85816/87353 Ethernet Camera (required for PEAK compatibility and remote support)
- Combine with other upgrades, installation
And onsite training

Onsite effort:

- Pre-configured to 'plug and play'
- Does not require SO Engineer on site
- Remote support included

How to contact us:

www.scientaomicron.com
info@scientaomicron.com