

Computer Ethernet Rack

Pre-configured and customized instrument PC With full support

Advantages:

- PEAK compatibility
- One year support including advice on plug-ins, e.g. manipulator control
- Win10-64 bit computer with hardware and software necessary for analyser control
- 'Plug and Play', customized and pre-installed
- Most recent SES software version with Real-time straight slit compensation
- Igor Pro 8 Software License



A Pre-configured and customized computer that meets the new PC requirements for PEAK compatibility. Latest platform; drivers and SES installed and tested. Including one year software support. Send instrument files to services-e-spec@Scientaomicron.com

Summary

Part numbers:

- 201042 Computer SES instrument ready
- 201043 Computer SES instrument ready DA30

Parts included:

- PC (does not include monitor/display)
- Instruction

Prerequisites:

- 102133 DOM upgrade
- 201040 Ethernet Rack upgrade (not required but recommended)
- 85827/85816 Ethernet Camera (not required but recommended)
- Depending on current configuration, additional upgrades may be required

Onsite effort:

- Pre-configured to 'plug and play'
- Does not require SO Engineer on site
- Remote support included

Options:

- Combine with other upgrades, installation
And onsite training

How to contact us:

www.scientaomicron.com
info@scientaomicron.com